





UP- AND DOWN- CONVERSION NANOPARTICLES


Up- and down-conversion nanoparticles are inorganic nanocrystals capable of converting near-infrared radiation to emission with large Stokes and/or anti-stokes shift. The upconversion photoluminescence mechanism includes the sequential absorption of several excitation photons through the long lifetime and ladder-like energy levels of trivalent lanthanide ions (ytterbium, erbium or thulium) embedded in an inorganic host matrix. The most efficient host matrix of UCNPs is NaYF₄ codoped with Yb³⁺, as a sensitizer, and Er³⁺ or Tm³⁺, as an activator. The luminescence of Lanthanide-doped nanoparticles is very stable without photobleaching and photoblinking issues. The overwhelming majority of the conversion nanoparticles prominent applications have originated from their unique luminescent properties. Today inorganic up- and down-conversion nanoparticles attract the interest in numerous fields, such as biomedical imaging, therapy agents, anti-counterfeit labels, 3D printing, information technology, solar cells, et al.



Anti-stokes fluorescence image upon excitation at 975 nm of TLS upconversion nanoparticles. A left cuvette filled with β -NaYF₄:Yb³⁺,Er³⁺/NaYF₄ and a right cuvette filled with β -NaYF₄:Yb³⁺,Tm³⁺/NaYF₄ demonstrate uniform green and UV/blue emission, respectively.





SEM image of upconversion nanoparticles with hexagonal prism shape.




TEM image of upconversion nanoparticle with hexagonal prism shape.

CONTACT US :

Te Linelo Systems BV
Mercurion 28 A
6903 PZ Zevenaar
The Netherlands
p. +31 316 340804
contact@tlsbv.nl



Stokes (left panel) and anti-Stokes (right panel) luminescent Spectrum of TLS upconversion nanoparticles with core/shell structure β -NaYF₄:Yb³⁺,Er³⁺/NaYF₄ irradiated with 975 nm laser light.



XRD pattern for TLS upconversion nanoparticles.

Product specification:

Composition	Up- and down- conversion nanoparticles		
Diameter	narrow dispersion	polydisperse	
	10-20 nm, 30-50 nm	10-100 nm	
Appearance	from clear to white solution		
Crystal formula	NaYREF ₄ , NaLaREF ₄ (RE: Yb, Er, Tm, Gd, Mn, Lu, Ce)		
Solvent	organic solvent or water		
Excitation wavelength	970-980 nm		
Sensitizer	Ytterbium (Yb ³⁺)		
Activator	Stokes emission	Anti-Stokes emission	
	Thulium (Tm ³⁺)	1700 nm	365/475/800 nm
	Erbium (Er ³⁺)	1530 nm	545/660 nm
Packaging	10 mg/mL in organic solvent in glass bottle, 1 – 10 mL		